


MANI[®] *VESSEL KNIFE*





MANI[®]

MANI® Vessel Knife (Blade & Handle) jointly invented by cardiovascular surgeon


Excellent sharpness

Samurai sword shape and cross-lapped edge based on fine processing technique provide ultimate sharpness.



Multiple combinations

Various angles can be made using two different angles of blades and handles.



Easiness and Simple Handling

Combination of compact edge and curved line blade provides elaborated scraping performance.



STERILE EO Sterilized using ethylene oxide

	Item code	Description	Package
	VKE100	Vessel Knife Blade 0 degree	5Pcs./Box
	VKE101	Vessel Knife Blade 35 degree	5Pcs./Box
	VKE500	Scraper Blade 0 degree	5Pcs./Box

Non-sterile

	Item code	Description	Package
	VKEH175-0	Vessel Knife Handle 0 degree	1Pc./Box
	VKEH175-10	Vessel Knife Handle 10 degree	1Pc./Box

MANI, INC.

8-3 Kiyohara Industrial Park, Utsunomiya, Tochigi, 321-3231, JAPAN
 Tel : +81-28-667-8346 Fax : +81-28-667-8305
 E-mail : surgical@ms.mani.co.jp URL : http://www.mani.co.jp