

MANI[®]

STAINLESS WIRE NEEDLE, STERILE



MANI[®]

MANI[®] Stainless Wire with

Y Cut Taper Needle


Better penetration, high ductility and tensile strength


Blunt Point Needle

High ductility, tensile strength and preventive measure against needle stick injury



■ STAINLESS WIRE NEEDLE, STERILE   0197

Needle Type: Y Point	USP Size	Description	Package	Item Code
	USP5	Wire Diameter 0.73~0.76mm Wire Length 450mm	2pcs/pack	0745SW-2
			3pcs/pack	0745SW-3
	USP6	Wire Diameter 0.83~0.86mm Wire Length 450mm	2pcs/pack	0845SW-2
			3pcs/pack	0845SW-3
	USP7	Wire Diameter 0.93~0.96mm Wire Length 450mm	2pcs/pack	0945SW-2
			3pcs/pack	0945SW-3

Needle Type: B Point	USP Size	Description	Package	Item Code
	USP5	Wire Diameter 0.73~0.76mm Wire Length 450mm	2pcs/pack	0745SWB-2
			3pcs/pack	0745SWB-3
	USP6	Wire Diameter 0.83~0.86mm Wire Length 450mm	2pcs/pack	0845SWB-2
			3pcs/pack	0845SWB-3
	USP7	Wire Diameter 0.93~0.96mm Wire Length 450mm	2pcs/pack	0945SWB-2
			3pcs/pack	0945SWB-3

Innovative Concept For Sternum Closure With *Blunt Point Needle*

**Penetrate
Through Bone!**



**Yet Safe
Against Glove!!**



■ Packaging 48pcs / Box



MANI, INC.

8-3 Kiyohara Industrial Park, Utsunomiya, Tochigi 321-3231 JAPAN
Tel : +81-28-667-8346 Fax : +81-28-667-8305
E-mail : surgical@ms.mani.co.jp URL : http://www.mani.co.jp