



MANI®



Disposable

MANI® ENT KNIFE

Disposable **MANI® ENT KNIFE**



*Example of use:
Posterior Nasal Neurectomy

- Small diameter shaft (2mm) is designed to secure the view and working space inside the nose.
- The angled blade tip provides better handling in creating incision in the nasal cavity.
- Mani® hard-fiber stainless steel creates ultimate sharpness and strength.
- Save your maintenance cost by using a disposable knife.



NKT35



Product Code	Product Name	Shaft	Packaging
NKT35	Nose Knife Triangle 3.5mm Straight	Straight	5Pcs/Box

NKT35C



Product Code	Product Name	Shaft	Packaging
NKT35C	Nose Knife Triangle 3.5mm Curved	Curved	5Pcs/Box

NKT35 NKT35C

*Actual size

STERILE EO Disposable

Do not use if package is damaged or opened.

Read instructions before use.

CE 0197

EC REP

Obelis s.a.
Bd. Général Wahis 53,
1030 Brussels, Belgium

Distributed by

Manufacturer :

MANI, INC.®

8-3 Kiyohara Industrial Park,
Utsunomiya, Tochigi, 321-3231 JAPAN
TEL: +81-28-667-7565 / FAX: +81-28-667-8305
URL: <http://www.mani.co.jp/en/>
E-mail: seg-ovs@ms.mani.co.jp

S01910(190712)