


MANI
STAINLESS WIRE NEEDLE,
STERILE



MANI

MANI Stainless Wire with

Y Cut Taper Needle


Better penetration, high ductility and tensile strength


Blunt Point Needle

High ductility, tensile strength and preventive measure against needle stick injury



STAINLESS WIRE NEEDLE, STERILE STERILEEO

Needle Type: Y Point	USP Size	Description	Package	Item Code
	USP5	Wire Diameter 0.73~0.76mm Wire Length 450mm	2pcs/pack	0745SW-2
			3pcs/pack	0745SW-3
	USP6	Wire Diameter 0.83~0.86mm Wire Length 450mm	2pcs/pack	0845SW-2
			3pcs/pack	0845SW-3
	USP7	Wire Diameter 0.93~0.96mm Wire Length 450mm	2pcs/pack	0945SW-2
			3pcs/pack	0945SW-3

Needle Type: B Point	USP Size	Description	Package	Item Code
	USP5	Wire Diameter 0.73~0.76mm Wire Length 450mm	2pcs/pack	0745SWB-2
			3pcs/pack	0745SWB-3
	USP6	Wire Diameter 0.83~0.86mm Wire Length 450mm	2pcs/pack	0845SWB-2
			3pcs/pack	0845SWB-3
	USP7	Wire Diameter 0.93~0.96mm Wire Length 450mm	2pcs/pack	0945SWB-2
			3pcs/pack	0945SWB-3

Innovative Concept For Sternum Closure With *Blunt Point Needle*

**Penetrate
Through Bone!**



**Yet Safe
Against Glove!!**



■ Packaging 48pcs / Box



Manufacturer :

MANI, INC.

8-3 KIYOHARA INDUSTRIAL PARK, UTSUNOMIYA, TOCHIGI, 321-3231, JAPAN

Phone: +81-28-667-7565 Fax: +81-28-667-8305

E-mail: seg-ovs@mani.inc URL: <http://www.mani.co.jp>