



MANI, INC.

8-3 KIYOHARA INDUSTRIAL PARK, UTSUNOMIYA, TOCHIGI, 321-3231, JAPAN

Exp.Sec./PHONE:81-28-667-2497 FAX:81-28-667-4964

ATT-24-0001

PAGE: 1 / 1

April 16th, 2024

Composition

Product – Suture Needle: MANI Eyeless Surgical Needle

We declare that the stainless-steel composition of the subject product is as follows.

The composition complies with JIS G 4308:2013 (equivalent to ASTM F899-20).

<MANI CODE X0>

[%]

C	Si	Mn	P	S	Ni	Cr
0.15 max	1.00 max	2.00 max	0.045 max	0.030 max	8.00~10.00	17.00~19.00

<MANI CODE X4>

[%]

C	Si	Mn	P	S	Ni	Cr
0.08 max	1.00 max	2.00 max	0.045 max	0.030 max	8.00~10.50	18.00~20.00