



JIZAI

— 自 在 —


020 04

Reduced risk of fracture
Reliable Assistance in root canals with severe calcifications, constriction, or curvature.


JIZAI has a high reputation for its centering ability and durability against fracture. Now, the newly added JIZAI 020 04 is designed to provide even a greater safety in treating difficult cases.

With the use of JIZAI 020 04, cutting load on JIZAI I (025 04) can be reduced by 56%!

● Cutting Image




● Cutting resistance test




Feature


JIZAI 020 04 with its "unique heat treatment and cross section", aims to minimize the risk of fracture and challenge complex cases.



Sequence



Specification

JIZAI 020 04	Size	Taper	Length	Quantity Per Unit
	020	04	21mm 25mm 28mm 31mm	3pcs

MANI, INC.

8-3 KIYOHARA INDUSTRIAL PARK, UTSUNOMIYA, TOCHIGI, 321-3231, JAPAN
Tel: +81-28-667-8592 E-Mail: dental.international@mani.inc
URL: www.mani.co.jp/en/